



Use Belden wire #8134 to connect PD Heads to the Box.

Color codes are as follows:

Color	PD Board	D-sub (male) pin#
Green	Test In	Unconnected
White	-15	3
Blue	+5	1
Red	+15	2
Brown	OUT	4
Black	GND	5
	PIN 9	9

Title ISS Remote PD Head		<i>LIGO Project</i> California Institute of Technology Massachusetts Institute of Technology		Cannot open file ligologo_1.jpg
Size: B	DCC Number: D030452-B-C			
Drawn by: Ben Abbott		Revision: B		
File: C:\Ben\Rickstuff\photodiode\PROTELSCH84 Date: 15/02/08 BUILT: 3:39:14 PM Sheet 1 of 1				