

D1900270-v7 LHAM7 WIP

6 April 2021

View is looking north/ -Y direction

12 vials are coming overnight to Matt from Liz. 4" x 4" areas will be used for FTIR sampling. Two areas will be FTIRed "before and after". That means the area will be FTIR before and then the area will be cleaned with high purity IPA (?) and then FTIRed again.

Assuming 12 vials: 1 control, 1 rinse and 10 vials. Map proposal

- 2-propyl + DI H<sub>2</sub>O*
- 1 & 2: Before and After floor, in the middle
  - 3 & 4: Before and After middle of north support tube surface (not in the holes)
  - 5: Swab around green area in support tube below nozzle
  6. Swab around convolutions of large bellows on the east door side
  7. Swab all support tube holes in the south support tube
  8. On the roof, to the west of the conflat
  9. Northeast corner of chamber
  10. Southeast corner of chamber (opposite to 9)

